

COPY OF MASTER FORMAT

CADILA HEALTHCARE LIMITED			Zydus Cadila
Annexure No.	0301-SOP-MFG-00666-01	Version No.	1.0, CURRENT
Annexure Title	Derivation of Force Control Limit in Compression Machine for "KN" Value		

Product Name					Batch No.	
Equipment ID.						
Average weight of tablets as per BMR (mg/g)	Avg. (mg/g)		Min. (mg/g)		Max. (mg/g)	

Parameter	Observed Limit	LHS	RHS
Compaction force at Target Average Weight	Observed Target Avg. weight (mg/g)		
	Observed Avg. compaction force – ‘KN’ (A)		
Compaction force at lower average weight	Observed Lower avg. weight (mg/g)		
	Observed Min. compaction force – ‘KN’ (B)		
	Observed Max. compaction force – ‘KN’ (C)		
Compaction force at higher average weight	Observed Higher average weight (mg/g)		
	Observed Min. compaction force – ‘KN’ (D)		
	Observed Max. compaction force – ‘KN’ (E)		

Parameter	Derived Limit	LHS	RHS
\$Derivation of main compaction force control limit (KN)	Lower limit: Average Compaction force at observed target Avg. weight (A) - up to 1 KN		
	Higher limit: Average Compaction force at observed target Avg. weight (A) + up to 1 KN		
@Derivation of rejection force control limit (KN)	Lower Limit (B): Min. value of individual compaction force observed at lower Avg. weight		
	Higher Limit (E): Max. value of individual compaction force observed at higher Avg. weight		
Done By:			
Verified By (Production):			
Verified By (IPQA): (Verify the set limit in SCADA against derived limits)			

Note: Ensure all the tablets comply to individual weight variation limit as per BMR specification.

\$ Compaction force can be changed by the operator during the compression run to maintain the desired target weight and other quality attributes of the product.

@ These are derived limit for each batch and fixed for complete batch run, which is set at initial batch setting stage.